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Powers KA, Kim PJ, Attinger CE, et al. Medstar Georgetown University Hospital, Center for Wound Healing.

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Early Experience With Negative Pressure Wound Therapy With Instillation in Acutely Infected Wounds

Kelly A Powers, DPM, MS Paul J Kim, DPM, MS Christopher E Attinger, MD John S Steinberg, DPM, Karen K Evans, MD
Zinnia M Rocha, BS Jesse R Smith, BA, MS Rex W Hung, BS

PURPOSE

This is a retrospective, historical cohort-controlled study assessing the effectiveness of Negative Pressure Wound Therapy (NPWT) versus Negative Pressure Wound Therapy with Instillation (NPWTi) in the treatment of infected wounds that require hospital admission and serial surgical debridement.

METHODS

Data obtained from electronic medical records from MedStar Georgetown University Hospital.

Study Dates:

- NPWT:
 - June, July, & August 2011
 - September, December 2011, & January 2012
- NPWTi:
 - June, July, & August 2012 (6 minute dwell time)
 - September, December 2012, & January 2013 (20 minutes dwell time)

Treatment Arms:

- Negative Pressure Wound Therapy (**NPWT**)
(Info V.A.C.® - KCI, San Antonio, TX)
- Negative Pressure Wound Therapy with Instillation (**NPWTi**)
(V.A.C. Ultra™ using Veraflo™ technology – KCI, San Antonio, TX)

Negative Pressure Setting:

- NPWT & NPWTi - 125mmHg Negative Pressure, Continuous

Instillation Solution:

- Prontosan™ (Polyhexanide + Betaine) (B Braun Medical Inc., Bethlehem, PA)
- 6 minute dwell (3.5 hours NPWT)
- 20 minute dwell time (2 hour NPWT)
- Volume varies. When foam is visibly saturated.

Subjects:


- June, July, & August 2011 (NPWT) N= 35
- June, July, & August 2012 (NPWTi – **6 Minute Dwell**) N= 34
- Sept, Dec 2011& Jan 2012 (NPWT) N= 39
- Sept, Dec 2012 & Jan 2013 (NPWTi) – **20 Minute Dwell**) N= 34

Variables:

- Number of O.R. Visits
- Length of Hospital Stay
- Time to Final Surgical Closure (*During the Admission Period*)
- Percent (%) Closed (*Prior to Discharge*)
- Percent (%) Improved - (Qualitative Cultures* - all subjects)
- Percent (%) Improved - (Qualitative Cultures* - excluding Gram negative bacteria, Corynebacterium, &/or Yeast)

* Qualitative Cultures based on OR visit #1 post-debridement qualitative cultures vs. OR visit #2 pre-debridement qualitative cultures.

INPATIENT INFECTED WOUND TREATMENT PROTOCOL


DEMOGRAPHICS


	NPWT	NPWTi
AGE	58 (26–83)	60 (20–88)
GENDER	M: 20 F: 15	M: 21 F: 18
RACE	A: 11 C: 20 H: 2 A: 1 Other: 1	A: 17 C: 16 H: 1 A: 0 Other: 5
BMI	42.2 (19.9–62.6)	31.5 (17.0–52.4)
Sept. & Dec. 2011 & Jan. 2012	63 (18–95)	55 (18–90)
Sept. & Dec. 2012 & Jan. 2013	60 (18–95)	55 (18–90)

WOUND LOCATION

NPWT

NPWTi

RESULTS

6 Minute Dwell Time

	NPWT N=35	NPWTi N=34	P
--	--------------	---------------	---

	NUMBER OF OR VISITS	2.9 ± 1.5	2.4 ± 0.9	*0.04
LENGTH OF STAY (days)	15.7 ± 10.8	11.9 ± 7.8	0.066	
TIME TO FINAL SURGICAL PROCEDURE (From Admission in days)	9.9 ± 6.3	7.8 ± 5.2	0.13	
PERCENT CLOSED (prior to discharge)	57% (20)	94% (32)	*.0001	
% IMPROVED: All Bacteria Included (Qualitative Cultures)	40% (14/35)	59% (20/34)	0.15	
% IMPROVED: Gram Negatives, Corynebacterium, Yeast Excluded (Qualitative Cultures)	64% (9/14)	90% (19/21)	0.09	

Improved= no growth or decrease in bacterial amount
Culture results were obtained from qualitative analysis
*p=0.05
Statistical significance calculation utilizing T-test of means, 2 tailed distribution.
Statistical significance calculation of % utilizing Fisher's exact test.

20 Minute Dwell Time

	NPWT N=39	NPWTi N=34	P
--	--------------	---------------	---

	NUMBER OF OR VISITS	3.0 ± 0.9	2.6 ± 0.9	*0.04
LENGTH OF STAY (days)	13.4 ± 4.8	11.4 ± 5.1	0.089	
TIME TO FINAL SURGICAL PROCEDURE (From Admission in days)	9.8 ± 4.3	7.5 ± 3.1	*0.012	
PERCENT CLOSED (prior to discharge)	67% (26)	80% (27)	0.29	
% IMPROVED: All Bacteria Included (Qualitative Cultures)	36% (14/39)	50% (17/34)	0.25	
% IMPROVED: Gram Negatives, Corynebacterium, Yeast Excluded (Qualitative Cultures)	62% (8/13)	65% (13/20)	1.0	

Improved= no growth or decrease in bacterial amount
Culture results were obtained from qualitative analysis
*p=0.05
Statistical significance calculation utilizing T-test of means, 2 tailed distribution.
Statistical significance calculation of % utilizing Fisher's exact test.

CONCLUSION

- ✓ Statistically significant decrease in the number of OR visits (6 & 20 Minute Dwell), % closed (6 Minute Dwell), time to final closure (20 Minute Dwell) using NPWTi as compared to NPWT.
- ✓ Trend for decrease in the length of hospital stay and higher % improvement in bacterial cultures using NPWTi as compared to NPWT.

DISCUSSION

- ✓ Negative pressure wound therapy with instillation shows promise for the adjunctive treatment of infected wounds that require hospital admission.
- ✓ Georgetown University Medical Center is leading a randomized, controlled prospective, multicenter study objectively examining the efficacy of NPWTi.

GEORGETOWN UNIVERSITY MEDICAL CENTER

Kelly A. Powers, DPM, MS – Research Fellow, Dept. of Plastic Surgery
Paul J. Kim, DPM, MS – Associate Professor, Dept. of Plastic Surgery
Christopher E. Attinger, MD – Professor, Dept. of Plastic Surgery
John S. Steinberg, DPM – Associate Professor, Dept. of Plastic Surgery
Karen K. Evans, MD – Assistant Professor, Dept. of Plastic Surgery
Zinnia M. Rocha, BS – Clinical Research Coordinator, Dept. of Plastic Surgery
Jesse R. Smith, BA, MS – Medical Student
Rex W. Hung, BS – Medical Student

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©2013 KCI Licensing, Inc. Todos los derechos reservados. Prontosan es una marca comercial de B. Braun Medical, Inc. Todas las demás marcas comerciales designadas aquí son propiedad de KCI Licensing, Inc., sus afiliadas y/o licenciatarios. DSL#13-0345.SAWC.P.LA • (Rev. 2/14)