Acelity™

EN

PREVENA™ INCISION MANAGEMENT SYSTEM

PREVENA™ PEEL & PLACE™ DRESSING

FOR USE WITH PREVENA™ THERAPY UNIT AND
KCI V.A.C.® THERAPY UNITS

APPLICATION INSTRUCTIONS FOR CLINICIANS ONLY

Rx Only


TABLE OF CONTENTS

Product Description and Indication for Use	5
Important Information for Users.....	5
Optimum Use Conditions.....	6
Contraindication.....	6
Warnings.....	6
Precautions	7
PREVENA™ Incision Management System Site Preparation.....	8
Drain Tubes and Pain Management Control Devices.....	8
PREVENA™ Incision Management System Dressings.....	9
PREVENA™ PEEL & PLACE™ System Kit Components	10
Dressing Application Instructions.....	11
Using the PREVENA™ Dressing with the PREVENA™ 125 Therapy Unit	12
Connecting the PREVENA™ Dressing to PREVENA™ 125 Therapy Units.....	12
Beginning Therapy.....	13
Duration of PREVENA™ Therapy.....	14
Dressing Removal.....	14
PREVENA™ 125 Therapy Unit Indicators and Alerts	15
PREVENA™ 125 Therapy Unit VISICHECK™ Feature	16
Correcting a Leak Condition.....	16
PREVENA™ 125 Therapy Unit Disposal.....	17
Battery Replacement on the PREVENA™ 125 Therapy Unit	18
Instructions For Patient	19
Daily Use	19
Sleeping	19
Showering and Bathing	19
Strenuous Activity.....	19
Cleaning.....	20
Device Disposal.....	20
Using the PREVENA™ Dressing with KCI V.A.C.® Therapy Units.....	20
Connecting the PREVENA™ Dressing to V.A.C.® Therapy Units	20

Setting Negative Pressure on the V.A.C. [®] Therapy Units	21
Alarm Resolutions	21
Specifications	21
Customer Contact Information	21
Electromagnetic Compatibility	22
Bibliography of Published Studies	26
Symbols Used	57

INSTRUCTIONS FOR USE

PREVENA™ INCISION MANAGEMENT SYSTEM

KCI customer contact information is located in the back of this guide.

PRODUCT DESCRIPTION AND INDICATION FOR USE

The PREVENA™ Incision Management System is intended to manage the environment of surgical incisions that continue to drain following sutured or stapled closure by maintaining a closed environment and removing exudate via the application of negative pressure wound therapy.

It consists of:

A PREVENA™ PEEL & PLACE™ Dressing and a source of negative pressure, which may be one of the following KCI therapy units:

- PREVENA™ 125 Therapy Unit
- ACTIV.A.C.™ Therapy Unit
- INFOV.A.C.™ Therapy Unit
- V.A.C. SIMPLICITY™ Therapy Unit
- V.A.C. FREEDOM™ Therapy Unit
- V.A.C.ULTA™ Therapy Unit
- V.A.C.RX4™ Therapy Unit

Clinical studies have been conducted on KCI Negative Pressure Incision Management Systems. Refer to the **Bibliography of Published Studies** in the back of this guide.

IMPORTANT INFORMATION FOR USERS

WARNING: The V.A.C.ULTA™, INFOV.A.C.™ and RX4™ Therapy Units are only indicated for use in acute care settings. Before transitioning the patient to home care, this therapy unit must be replaced with one indicated for home use, such as the PREVENA™ 125, ACTIV.A.C.™, V.A.C. FREEDOM™ or V.A.C. SIMPLICITY™ Therapy Unit.

CAUTION: The PREVENA™ Therapy System should be applied and removed only by qualified physicians or nurses.

As with any prescription medical device, failure to carefully read and follow all instructions and safety information prior to use may lead to improper product performance.

The PREVENA™ Therapy System dressings and therapy unit canisters are disposable and are for single use only. Re-use of disposable components may result in wound contamination, infection and/or failure of the wound to heal.

WARNING: DO NOT use with V.A.C. VERAFL™ Therapy (Instillation) provided by the V.A.C.ULTA™ Therapy Unit. Instillation into the incision site may result in pooling of fluid which may result in maceration.

For pressure settings and connection information for use of the PREVENA™ Dressings with the V.A.C.® Therapy Units listed in the **Product Description and Indication for Use section** above, see the **Using the PREVENA™ Dressing with KCI V.A.C.® Therapy Units** section.

OPTIMUM USE CONDITIONS

For maximum benefit the PREVENA™ Therapy System should be applied immediately post surgery to clean surgically closed wounds. It is to be continuously applied for a minimum of two days up to a maximum of seven days. It can transition home with the patient.

The PREVENA™ Therapy System will not be effective in addressing complications associated with the following:

- ischemia to the incision or incision area
- untreated or inadequately treated infection
- inadequate hemostasis of the incision
- cellulitis of the incision area

The PREVENA™ Therapy System should not be used to treat open or dehisced surgical wounds. The V.A.C.® Therapy System should be considered for treatment of these wounds.

A V.A.C.® Therapy Unit with the smallest available canister should be considered for patients who have exudate volumes that exceed the PREVENA™ 45ml Canister.

The PREVENA™ Therapy System should be used with caution in the following patients:

- patients with fragile skin surrounding the incision as they may experience skin or tissue damage upon removal of the PREVENA™ Dressing
- patients who are at an increased risk of bleeding from the incision associated with the use of anticoagulants and/or platelet aggregation inhibitors

CONTRAINDICATION

- sensitivity to silver

WARNINGS

The PREVENA™ Therapy System is not intended to manage open or dehisced wounds.

DO NOT use with V.A.C. VERAFL™ Therapy (Instillation) provided by the V.A.C.ULTA™ Therapy Unit. Instillation into the incision site may result in pooling of fluid which may result in maceration.

Bleeding: Before applying the PREVENA™ Therapy System to patients who are at risk of bleeding complications due to the operative procedure or concomitant therapies and/or co-morbidities, ensure that hemostasis has been achieved and all tissue planes have been approximated. If active bleeding develops suddenly or in large amounts during therapy, or if frank blood is seen in the tubing or in the canister, the patient should leave the PREVENA™ Dressing in place, turn off the therapy unit and seek immediate emergency medical assistance.

Infected Wounds: As with any wound treatment, clinicians and patients/caregivers should frequently monitor the patient's wound, periwound tissue and exudate for signs of infection or other complications. Some signs of infection are fever, tenderness, redness, swelling, itching, rash, increased warmth in the wound or periwound area, purulent discharge or strong odor. Infection can be serious, and can lead to complications such as pain, discomfort, fever, gangrene, toxic shock, septic shock and/or fatal injury. Some signs or complications of systemic infection are nausea, vomiting, diarrhea, headache, dizziness, fainting, sore throat with swelling of the mucus membranes, disorientation, high fever, refractory and/or orthostatic hypotension or erythroderma (a sunburn-like rash). Silver in the interface layer of the PREVENA™ Dressing is not intended to treat infection, but to reduce bacterial colonization in the fabric. **If infection develops, the PREVENA™ Therapy System should be discontinued until the infection is treated.**

Allergic Response: The PREVENA™ Dressing has an acrylic adhesive coating and a skin interface layer with silver, which may present a risk of an adverse reaction in patients who are allergic or hypersensitive to acrylic adhesives or silver. If a patient has a known allergy or hypersensitivity to these materials, do not use the PREVENA™ Therapy System. If any signs of allergic reaction, irritation or hypersensitivity develop, such as redness, swelling, rash, urticaria or significant pruritus, patient should consult a physician immediately. If bronchospasm or more serious signs of allergic reaction appear, the patient should turn off the therapy unit and seek immediate emergency medical assistance.

Defibrillation: Remove the PREVENA™ Dressing if defibrillation is required in the area of dressing placement. Failure to remove the dressing may inhibit transmission of electrical energy and/or patient resuscitation.

Magnetic Resonance Imaging (MRI): All KCI therapy units, including the PREVENA™ 125 Therapy Unit, are MR unsafe. Do not take therapy units into the MR environment. The PREVENA™ Dressings can typically remain on the patient with minimal risk in an MR environment. Interruption of PREVENA™ Therapy during MRI may reduce the effectiveness of PREVENA™ Therapy. The PREVENA™ Therapy System Dressings pose no known hazards in an MR environment with the following conditions of use: static magnetic field of 3 Tesla or less, spatial gradient field of 720 Gauss/cm or less and maximum whole-body-averaged specific absorption rate (SAR) of 3W/kg for 15 minutes of scanning.

Diagnostic Imaging: The PREVENA™ Dressing contains metallic silver that may impair visualization with certain imaging modalities.

Hyperbaric Oxygen Therapy (HBO): Do not take therapy units or PREVENA™ Dressings into a hyperbaric oxygen chamber. They are not designed for this environment and **should be considered a fire hazard**. If the PREVENA™ Therapy System is reinitiated after HBO treatment, do not readhere the same dressing; a new dressing must be applied.

Canister Full: If at any time while using the PREVENA™ Therapy System the canister becomes full of fluid other than blood, indicated by a therapy unit alert or visual inspection, the patient should turn off the therapy unit and contact the treating physician.

Standard Operation: Do not use accessories or materials not provided with the PREVENA™ Therapy System. For a list of acceptable therapy units with which PREVENA™ Dressings may be used, see the **Product Description and Indication for Use** section.

PRECAUTIONS

Standard Precautions: To reduce the risk of transmission of bloodborne pathogens, apply standard precautions for infection control with all patients, per institutional protocol, regardless of their diagnosis or presumed infection status.

Circumferential Dressing Application: Avoid applying the PREVENA™ Dressing circumferentially. In cases where the clinician determines that the benefits of applying the PREVENA™ Dressing circumferentially outweigh the risk of circulatory compromise, extreme care should be taken not to stretch or pull the dressing when securing it. Attach the dressing loosely and stabilize edges with an elastic wrap if necessary. It is crucial to systematically and recurrently palpate distal pulses and assess distal circulatory status. If circulatory compromise is suspected, discontinue therapy and remove dressing.

Electrodes or Conductive Gel: Do not allow the PREVENA™ Dressing to come in contact with EKG or other electrodes or conductive gels during electronic monitoring or when taking electronic measurements.

Dressing Components: The PREVENA™ Dressing contains ionic silver (0.019%). Application of products containing silver may cause temporary tissue discoloration.

- Always use PREVENA™ Dressings and canisters from sterile packages that have not been opened or damaged.
- All dressing components and canisters of the PREVENA™ Therapy System are for single use only. Do not re-use any component of this system.
- To avoid trauma to the skin, do not pull or stretch the adhesive border of the dressing during application.

Compressive Garments or Dressings: Avoid tight compressive garments or dressings (such as surgical bras, elastic bandage wraps or abdominal binders) to prevent forcibly pressing the PREVENA™ Dressing into soft tissue.

PREVENA™ INCISION MANAGEMENT SYSTEM SITE PREPARATION

1. Prior to surgery, shave or clip the surgical area where the dressing will be applied to improve dressing adhesion and seal integrity.
2. Gather all items needed for application:
 - sterile wound cleaning solution, e.g. water, saline or alcohol
 - sterile gauze or other material to clean application site
 - all components of the PREVENA™ Therapy System (dressing and therapy unit)
3. After surgery, cleanse the application site with sterile gauze and sterile wound cleaning solution using a circular motion beginning at the center of the surgical area and extending outward to ensure that the site is free of foreign material.
4. Pat the application site dry with sterile gauze. To ensure proper adhesion, the application site must be completely dry before dressing is applied.

DRAIN TUBES AND PAIN MANAGEMENT CONTROL DEVICES

The PREVENA™ Therapy System can be used with both drain tubes and pain devices, provided the dressing is not placed over tubing where it exits the skin. Surgical drains must be routed under the skin beyond the boundary of the dressing and function independently of the PREVENA™ Therapy System.

NOTE: While the concomitant use of surgical drains is allowable with the PREVENA™ Therapy System, the system must not be used as an outlet or reservoir for the drain.

PREVENA™ INCISION MANAGEMENT SYSTEM DRESSINGS


The PREVENA™ Therapy System features a PREVENA™ PEEL & PLACE™ Dressing. This dressing contains a skin interface layer that includes 0.019% ionic silver. Silver in the interface layer of the PREVENA™ Dressing is not intended to treat infection, but to reduce bacterial colonization in the fabric. *In vitro* log reduction tests, conducted without application of negative pressure, exposed samples of the skin interface layer to a six log challenge of each of the microorganisms listed below. Following inoculation, samples were tested for microbial counts immediately (day 0) and after incubation at 90°F (32°C) in diluted nutrient broth for 1, 3, 5 and 7 days.

The log reductions from the day 0 values are provided in the table below.


Challenge Organism	Mean Log Reduction from Day 0			
	Day 1	Day 3	Day 5	Day 7
<i>Escherichia coli</i> (ATCC 8739)	2.2	4.0	3.9	4.5
<i>Pseudomonas aeruginosa</i> (ATCC 09027)	2.0	3.9	3.5	3.7
<i>Staphylococcus aureus</i> (ATCC 6538)	1.6	3.6	3.6	3.5
<i>Klebsiella pneumonia</i> (ATCC 4352)	1.4	1.8	2.7	3.5
<i>Candida albicans</i> (ATCC 10231)	2.5	3.1	3.2	3.2
<i>Aspergillus niger</i> (ATCC 16404)	2.2	4.1	4.0	3.6

PREVENA™ PEEL & PLACE™ SYSTEM KIT COMPONENTS


The PREVENA™ PEEL & PLACE™ System Kit contains the following single use, disposable components.

PREVENA™ PEEL & PLACE™ Dressing - 13cm - a specially designed dressing for application to the surgical area (**13 cm kit only**)


PREVENA™ PEEL & PLACE™ Dressing - 20cm - a specially designed dressing for application to the surgical area (**20 cm kit only**)


PREVENA™ Patch Strips - used to help seal leaks around dressing

PREVENA™ Therapy V.A.C.® Connector - used to connect the PREVENA™ Dressing to an approved KCI V.A.C.® Therapy Unit

Ruler - the removable label may be used as needed to record the date of dressing application or removal

PREVENA™ 45ml Canister - a sterile reservoir for collection of wound fluids

PREVENA™ 125 Therapy Unit - delivers negative pressure to the surgical area. The unit is battery powered. The non-sterile PREVENA™ 125 Therapy Unit Carrying Case is provided to facilitate patient mobility.


DRESSING APPLICATION INSTRUCTIONS

(Illustrations in the steps below show PREVENA™ PEEL & PLACE™ Dressing - 20cm)


CAUTION: If the dressing covers the umbilicus, the umbilicus must first be fully filled with an anti-microbial petroleum gauze prior to dressing application.

WARNING: DO NOT use with V.A.C. VERAFL™ Therapy (Instillation) provided by the V.A.C. ULTA™ Therapy Unit. Instillation into the incision site may result in pooling of fluid which may result in maceration.

1. Open the sterile dressing package and remove dressing and patch strips using aseptic technique. Do not use if package has been torn or the sterile seal has been compromised.


2. Gently peel back the center strip on the back of the dressing exposing the pull tabs and adhesive.

3. Center and apply the dressing over the closed wound or incision ensuring that the adhesive will not contact or cover the surgical closure. Orient the dressing on the patient to eliminate sharp bends or kinks in the tubing.

4. Remove the remaining bottom adhesive covers by grasping the bottom tabs and gently pulling.

5. Firmly press around the dressing to ensure a good seal where the adhesive contacts the skin.


6. Remove top stabilization layers.

7. Connect to therapy unit.

For connection to the PREVENA™ 125 Therapy Unit, see the **Connecting the PREVENA™ Dressing to PREVENA™ 125 Therapy Units** section.

For connection to other KCI V.A.C.® Therapy Units, see the **Connecting the PREVENA™ Dressing to V.A.C.® Therapy Units** section.

For Dressing Removal, see the Dressing Removal section.


NOTE: The removable label on the supplied ruler may be used as needed to record date of dressing application or removal.


USING THE PREVENA™ DRESSING WITH THE PREVENA™ 125 THERAPY UNIT

CONNECTING THE PREVENA™ DRESSING TO PREVENA™ 125 THERAPY UNITS

1. Remove the PREVENA™ 45ml Canister from the sterile package. Do not use if package has been torn or the sterile seal has been compromised.


2. Insert the canister into the PREVENA™ 125 Therapy Unit and slide inward until canister clicks. Canister is fully inserted when the side tabs are flush with the body of the therapy unit.

3. Connect the dressing tubing to the canister tubing by twisting the connectors until they lock.

4. Begin therapy.

BEGINNING THERAPY


1. Press and hold the On/Off button for two seconds; an audible beep will confirm that therapy is on. A green light on the front of the unit indicates that therapy is on.


NOTE: Pressing the On/Off button begins the 192 hour (eight day expected service life) life cycle of the therapy unit. Turning the therapy unit off stops the life cycle counter. Turning the therapy unit on for purposes other than delivering therapy reduces the life cycle of the therapy unit. It is not recommended to press the On/Off button until therapy is ready to begin.

NOTE: To turn therapy unit off, press and hold the On/Off button for five seconds.


2. With therapy on, assess dressing to ensure integrity of seal.


- ✓ • The dressing should have a wrinkled appearance and the foam bolster should be compressed.

- ✗ • If the foam bolster is not compressed or the therapy unit alerts, see the **PREVENA™ 125 Therapy Unit Indicators and Alerts** section.

3. Place the therapy unit into the PREVENA™ 125 Therapy Unit Carrying Case. Make sure the display is visible through the opening in the carrying case when the front flap is lifted.

4. The PREVENA™ 125 Therapy Unit Carrying Case has an integrated belt loop and a separate adjustable strap to allow for versatile positioning.

CAUTION: Do not wrap carrying case strap or dressing tubing around neck.



DURATION OF PREVENA™ THERAPY

- Therapy should be continuous for a minimum of two days up to a maximum of seven days.
- The PREVENA™ 125 Therapy Unit will automatically time-out after 192 hours (eight days) of cumulative run time.
- Patients should be instructed to contact their treating physician and not to turn therapy off unless:
 - advised by the treating physician
 - bleeding develops suddenly or in large amounts during therapy
 - there are signs of allergic reaction or infection
 - the canister is full of fluid
 - batteries need to be changed
 - system alerts must be addressed
- Patient should be instructed to contact the treating physician if therapy unit turns off and cannot be restarted before therapy is scheduled to end, or if canister becomes full of fluid.
- At end of therapy, patient should return to treating physician for dressing removal.

DRESSING REMOVAL

NOTE: If dressing is lifted to observe wound, do not re-adhere the same dressing; a new dressing must be applied.

WARNING: Dressings should always be removed in-line with the sutures and NEVER across the sutures.

1. Turn the therapy unit off by pressing and holding the On/Off button.
2. Gently stretch the drape/dressing horizontally to release the adhesive from the skin. Do not peel vertically. Remove the drape/dressing in-line with the sutures, NEVER across the sutures.
3. Clean any residual adhesive.

If a new dressing is to be applied:

1. Ensure the incision area is clean by using an alcohol swab or an antiseptic wipe.
2. Allow skin to dry completely before applying a new drape.
3. Follow the **PREVENA™ PEEL & PLACE™ Dressing Application** section.

PREVENA™ 125 THERAPY UNIT INDICATORS AND ALERTS

Caution Indicator

Visual Alerts - Solid LEDs cannot be turned off by the user. Visual alerts will only stop when the alert condition has been corrected.

Audible Alarm (Muted)

Audible Alerts - Repeated beeps (which in some cases will increase in volume) can be temporarily muted (paused) by pressing the On/Off button once. The audible alert will re-occur unless the alert condition has been corrected.

The PREVENA™ 125 Therapy Unit will provide audible and visual alerts as shown below.

Battery Level

If alert conditions cannot be corrected, patient should contact treating physician. For further product support see **Customer Contact Information** section.

Leak Alert

Audible Alert can be muted for eight hours.

One beep, one solid yellow light.

See **Correcting a Leak Condition** section.

When a leak condition is corrected the therapy unit will cancel alert. There may be a delay between when the leak is corrected and when the alert is discontinued.

Canister Full Alert

Audible Alert can be muted for one hour.

Two beeps, one solid yellow light.

Visually inspect canister. If full or near full, turn therapy unit off and call the treating physician immediately.

Low Battery Alerts

Audible Alert can be muted for one hour.

LOW - One slow beep, one solid yellow light. Be prepared to change batteries.

CRITICAL - One beep, repeating rapidly, increasing in volume, one solid yellow light. Change batteries immediately.

Battery replacement will cancel the alert. See **Battery Replacement on the PREVENA™ 125 Therapy Unit** section.

System Error Alert

Audible Alert can be muted for one hour.

One beep, repeating rapidly, increasing in volume, two solid yellow lights. Turn the unit off and then on again.

If alert continues, notify the treating physician.

Device Lifecycle Expired

Three solid yellow lights. If the lifecycle expires while therapy unit is on, a beep will sound for 15 seconds and then automatically shut off.

If the therapy unit is off, has exceeded the run time, and an attempt is made to turn the therapy unit on, the therapy unit will alert for three seconds and automatically shut off. Call the treating physician.

PREVENA™ 125 THERAPY UNIT VISICHECK™ FEATURE

Best

Good

Marginal

To ensure proper application of the PREVENA™ Dressing, the PREVENA™ 125 Therapy Unit provides the VISICHECK™ Feature.

By double pressing the On/Off button, the unit will display the leak rate of the system for three seconds.


To prevent nuisance leak alarms, the leak rate status should be "Best" (one light illuminated) or "Good" (two lights illuminated).

When the therapy unit detects a leak:


- A visual and audible leak alert will activate. See the **PREVENA™ 125 Therapy Unit Indicators and Alerts** section.
- Therapy unit will turn on more frequently.

CORRECTING A LEAK CONDITION

(Illustrations in the steps below show PREVENA™ PEEL & PLACE™ Dressing - 20cm)


1. With therapy unit on, slowly press firmly around dressing edge to ensure good contact between adhesive and skin.


2. If a leak is identified, use PREVENA™ Patch Strips (located in dressing package) to help seal leaks around dressing. If large wrinkles are present, place patch strips so they run in line along the length of the wrinkle and not across the wrinkle.

Check Canister and Tubing Connection


1. Ensure canister is securely locked in the therapy unit. When canister is installed, a distinct click will be heard indicating it has been properly installed. Side tabs on canister should be flush with unit.


2. Check tubing connectors to ensure they are fully engaged and locked.

Indications That a Leak Condition Has Been Corrected

- Therapy unit will become quiet, running only intermittently.
- Audible leak alert will stop; visual alert will turn off. There may be a brief delay between when the leak is corrected and the alert is discontinued.
- The PREVENA™ Dressing will be compressed.

Dressing compressed - system pressure acceptable.

Dressing not compressed - system pressure not acceptable. Return to the **Correcting a Leak Condition** section to continue pressure correction steps.

PREVENA™ 125 THERAPY UNIT DISPOSAL



At the end of therapy, the patient should return the PREVENA™ 125 Therapy Unit to the physician for disposal. Dispose of all waste according to local requirements. Improper disposal may run the risk of regulatory non-compliance.

BATTERY REPLACEMENT ON THE PREVENA™ 125 THERAPY UNIT

Battery replacement should be done as quickly as possible after a low battery alert to prevent therapy down time.


1. Turn therapy unit off (press and hold On/Off button for five seconds).
2. Remove therapy unit from carrying case. Turn the therapy unit over to expose the back side of the unit.

3. Locate the battery door and push to slide open. Install three AA batteries (lithium batteries are recommended for optimal performance) into the battery compartment.

NOTE: Always replace with new batteries. Do not mix new batteries with used batteries.

NOTE: The inside of the battery compartment is stenciled with "+" positive and "-" negative to aid in proper battery installation.

4. Close battery door.
5. Return therapy unit to carrying case.
6. Turn therapy unit on to resume therapy (press and hold On/Off button for two seconds).

INSTRUCTIONS FOR PATIENT

Review the following information with the patient prior to discharge. This information is summarized in the PREVENA™ Therapy System Patient Guide which must be provided to the patient at discharge.

WARNING: Inform the patient that the V.A.C.ULTA™, INFOV.A.C.™ and RX4™ Therapy Units should not be sent home. They should be sent home with another KCI negative pressure therapy unit such as the PREVENA™ 125, ACTIV.A.C.™, V.A.C. FREEDOM™ or V.A.C. SIMPLICITY™ Therapy Unit.

DAILY USE

The PREVENA™ 125, ACTIV.A.C.™, V.A.C. FREEDOM™ and V.A.C. SIMPLICITY™ Therapy Units are portable and small enough that they may be worn beneath clothing during normal patient activities as approved by the treating physician.

CAUTION: Advise patient to NOT SUBMERGE therapy unit or dressing in liquid and to ensure therapy unit is not pulled into a tub or sink where it may become submerged.

CAUTION: The PREVENA™ 125 Therapy Unit is a medical device not a toy. Keep away from children, pets and pests as they can damage the dressing and therapy unit and affect performance. Keep therapy unit free of dust and lint.

SLEEPING

Instruct patient to:

- place the therapy unit in a position where tubing will not become kinked or pinched.
- ensure therapy unit will not be pulled off a table or fall to the floor during sleep.

SHOWERING AND BATHING

Advise patient of the following recommendations:

- Light showering is permissible, bathing is not.
- While showering, the device and dressing should be protected from prolonged direct spray and/or being submerged.
 - The PREVENA™ 125 Therapy Unit may be hung from a soap/shampoo holder or shower head if protected from prolonged direct spray.
 - If using a KCI V.A.C.® Therapy Unit, disconnect the unit from the dressing before bathing or showering. Refer to therapy unit's user manual for instructions.
 - Dressing may be exposed to common shower soaps and rinsed with indirect shower stream. Do not submerge dressing. Do not remove dressing.
- When towel drying, avoid disturbing or damaging dressing.

STRENUOUS ACTIVITY

Advise patient as to when and at what level physical activities may be resumed. It is recommended that patients avoid strenuous activity while using the therapy system.

CLEANING

Advise patient that the therapy unit and therapy unit carrying case can be wiped with a damp cloth using a mild household soap solution that does not contain bleach.

DEVICE DISPOSAL

At the end of therapy, the patient should return therapy unit to physician for disposal. Dispose of all waste according to local requirements. Improper disposal may run the risk of regulatory non-compliance.

USING THE PREVENA™ DRESSING WITH KCI V.A.C.® THERAPY UNITS

When directed by the treating physician, PREVENA™ Dressings can be used with negative pressure wound therapy provided by INFOV.A.C.™, RX4™, V.A.C. SIMPLICITY™, ACTIV.A.C.™, V.A.C. FREEDOM™, and V.A.C.ULTA™ Therapy Units.

WARNING: Refer to the therapy unit's user manual for complete instructions for use and safety information before initiating therapy.

WARNING: DO NOT use with V.A.C. VERAFL™ Therapy (Instillation) provided by the V.A.C.ULTA™ Therapy Unit. Instillation into the incision site may result in pooling of fluid which may result in maceration.

WARNING: The V.A.C.ULTA™, INFOV.A.C.™ and RX4™ Therapy Units are only indicated for use in acute care settings. Before transitioning the patient to home care, this therapy unit must be replaced with one indicated for home use, such as the PREVENA™ 125, ACTIV.A.C.™, V.A.C. FREEDOM™ or V.A.C. SIMPLICITY™ Therapy Unit.

NOTE: When PREVENA™ Dressings are used with V.A.C.® Therapy, the SENSAT.R.A.C.™ wound pressure sensing function on the therapy unit will be disabled. Instead, the PREVENA™ Dressing foam bolster compression state should be used to verify that incision site pressure is acceptable.

NOTE: Consider using the smallest available canister for the selected V.A.C.® Therapy Unit.

CONNECTING THE PREVENA™ DRESSING TO V.A.C.® THERAPY UNITS

A KCI PREVENA™ Therapy V.A.C.® Connector will be needed to connect the PREVENA™ Dressing to a KCI V.A.C.® Therapy Unit. This connector, which is available in the PREVENA™ Dressing package, **must** be used for negative pressure wound therapy to work effectively and accurately. No other method for connecting the dressing to the therapy unit should be attempted.

1. Connect the PREVENA™ Dressing tubing to the V.A.C.® Therapy Unit canister tubing using the PREVENA™ Therapy V.A.C.® Connector:
 - Push the connectors together.
 - Twist the connectors to lock.

2. Ensure clamp on canister tubing is open.

SETTING NEGATIVE PRESSURE ON THE V.A.C.® THERAPY UNITS

Activate V.A.C.® Therapy at -125mmHg continuous. Do not choose any other negative pressure setting, intermittent or DPC modes of negative pressure.

For V.A.C.ULTA™ Therapy Units:

- Select V.A.C.® Therapy **only**.
- Do NOT select VERAFLÓ™ Therapy (see the previous **WARNING**).
- See the **Choose Therapy** section of the V.A.C.ULTA™ User Manual for more information.

ALARM RESOLUTIONS

KCI V.A.C.® Therapy Unit alarms should be addressed in a timely manner. Refer to the appropriate therapy unit user manual for complete information on alarm resolutions. Refer to the **Correcting a Leak Condition** section for correcting a leak in the dressing.

SPECIFICATIONS

Storage Conditions

Temperature Range - 4°F (-20°C) to 140°F (60°C)
Relative Humidity Range 15% - 95%, non-condensing

Operating Conditions

Temperature Range 41°F (5°C) to 104°F (40°C)
Atmospheric Pressure Range 1060 hpa (-1253 ft / -381.9 m)
(for Optimum Performance) to 700 hpa (9878 ft / 3010 m)

Equipment not suitable for use in presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.

The PREVENA™ Therapy System is classified as a Type BF applied part under IEC 60601-1.

IP24 - Protection against solid objects greater than 12.5 mm and against liquid water sprays for short periods of time.

All alerts are classified as low priority according to IEC 60601-1-8.

Conforms to IEC 60601-1, IEC 60601-1-6, IEC 62366, IEC 60601-1-8, IEC 60601-1-11.

CUSTOMER CONTACT INFORMATION

For questions regarding this product, supplies, maintenance or additional information about KCI products and services, please contact KCI or a KCI authorized representative, or:

In the US call 1-800-275-4524 or visit www.acelity.com

KCI USA, Inc. 12930 IH10 West, San Antonio, TX 78249

ELECTROMAGNETIC COMPATIBILITY

All electrical equipment may produce interference. If interference is suspected, move equipment away from sensitive devices or contact the manufacturer.

The PREVENA™ 125 Therapy Unit needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information in the following tables.

Guidance and manufacturer's declaration - electromagnetic emissions			
The PREVENA™ 125 Therapy Unit is intended for use in the electromagnetic environment specified below. The customer or user of the PREVENA™ 125 Therapy Unit should assure that it is used in such an environment.			
Emission Test	Compliance	Electromagnetic environment	
RF emissions CISPR 11	Group 1	The PREVENA™ 125 Therapy Unit uses RF energy only for internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR 11	Class B	The PREVENA™ 125 Therapy Unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic emissions IEC 61000-3-2	Not applicable	Battery operated device	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable		
WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.			
Portable and mobile RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the PREVENA™ 125 Therapy Unit. Otherwise, degradation of the performance of this equipment could result. More precisely, the minimum recommended separation distance should be calculated from the equation applicable to the frequency of the transmitter, as noted in the guidance below.			
NOTE: This equipment has been tested and found to comply with the limits for medical devices to IEC 60601-1-2: 2014 4th edition. These limits and test levels are intended to provide reasonable safety with regard to electromagnetic disturbances when the device is used in a typical medical installation.			

Guidance and manufacturer's declaration - electromagnetic immunity

The PREVENA™ 125 Therapy Unit is intended for use in the electromagnetic environment specified below. The customer or user of the PREVENA™ 125 Therapy Unit should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not Applicable	Battery operated device
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not Applicable	Battery operated device
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_r (>95 % dip in U_r) for 0.5 cycle 40 % U_r (60 % dip in U_r) for 5 cycles 70 % U_r (30 % dip in U_r) for 25 cycles <5 % U_r (>95 % dip in U_r) for 5 cycles	Not Applicable	Battery operated device
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m 50Hz or 60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_r is the a.c. mains voltage prior to application of the test level.

Recommended separation distances between portable and mobile RF communications equipment and the PREVENA™ 125 Therapy Unit

The PREVENA™ 125 Therapy Unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the PREVENA™ 125 Therapy Unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the PREVENA™ 125 Therapy Unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter in watts (W)	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 0.35\sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7\sqrt{P}$
0.01	.12	.04	.07
0.1	.38	.11	.22
1	1.2	.35	.7
10	3.8	1.1	2.2
100	12	3.5	7

For transmitters rated at a maximum output power not listed above, the recommended separate distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1, At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from surfaces, objects and people.

Guidance and manufacturer's declaration - electromagnetic immunity

The PREVENA™ 125 Therapy Unit is intended for use in an electromagnetic environment specified below. The customer or user of the PREVENA™ 125 Therapy Unit should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not Applicable	Portable and mobile RF communications equipment should be used no closer to any part of the PREVENA™ 125 Therapy Unit than the recommended separation distance calculated from the equation application to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz	<p>Recommended separation distance</p> $d = 0.35\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 0.7\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m)</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey¹, should be less than the compliance level in each frequency range².</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE 1, At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

¹ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the PREVENA™ 125 Therapy Unit is used exceeds the applicable RF compliance level above, the PREVENA™ 125 Therapy Unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the PREVENA™ 125 Therapy Unit.

² Over the frequency range 150kHz, field strengths should be less than 3 V/m.

BIBLIOGRAPHY OF PUBLISHED STUDIES

- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" *Surgical Innovations*. 2009 June 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." *International Wound Journal*. 2011 Feb;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." *Techniques in Foot & Ankle Surgery*. 2013 December 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." *American Journal of Surgery*. 2013 January 30.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." *International Wound Journal*. 2013 Feb 19. [Epub ahead of print]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." *Diseases of the Colon & Rectum*. 2013 December 1;56(12):1403-8
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." *Journal of Cardiothoracic Surgery*. 2011 December 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." *Annals of Plastic Surgery* 2013 Oct 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." *Journal of Wound Care*. 2012 Aug 21(8):386-8.
- Faroqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." *Journal of Diabetic Foot Complications*. 2013 September 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." *Wounds*. 2012 Nov 24(11):308-11
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." *Journal of Orthopaedic Trauma*. 2006 Nov-Dec;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." *Journal of Thoracic and Cardiovascular Surgery*. 2013 May;145(5):1387-92.
- Haghshenasskashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." *The British Journal of Diabetes and Vascular Disease*. 2011 Jan/Feb;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." *Eplasty* 2013 September 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." *Bone and Joint Research*. 2013 December 1;2(12):276-84.
- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." *Wound Repair and Regeneration*. 2011 Sep;19(5):588-96.

- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." *Surgical Innovations*. 2013 September 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." *Annals of Surgery*. 2012 Jun;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." *Journal of Vascular Surgery*. 2013 Mar;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." *International Orthopaedics*. 2012 Apr;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." *Surgical Infections*. 2013 Jun;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." *Journal of Surgical Orthopaedic Advances*. 2010 Jun;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." *Journal of Trauma*. 2006 Jun;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." *Journal of Orthopedic Trauma*. 2012 Jan;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." *Journal of Plastic, Reconstructive, and Aesthetic Surgery*. 2013 Mar;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." *American Journal of Surgery*. 2012 Dec;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." *Surgical Innovations*. 2012 Mar;19(1):67-75.

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-
- Two small diagrams illustrating the connection process. The left diagram shows two connectors being joined together with arrows indicating the direction of assembly. The right diagram shows a single connector attached to a tube, with an arrow pointing to the pinza (clip) mechanism on the side of the connector.
1. Conecte el tubo del apósito PREVENA™ para la terapia V.A.C.® al tubo del contenedor de la unidad de terapia V.A.C.® mediante el conector de V.A.C.® de la terapia PREVENA™.
 - Presione los conectores entre sí.
 - Gire los conectores para bloquearlos.
 2. Asegúrese de que la pinza en el tubo del contenedor esté abierta.

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- Atkins BZ y cols. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" *Surgical Innovations*. Junio de 2009; 16(2):140-6.
- Atkins BZ y cols. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy". *International Wound Journal*. Febrero de 2011; 8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures". *Techniques in Foot & Ankle Surgery*. 1 de diciembre de 2013;12(4):172-4.
- Blackham AU y cols. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy". *American Journal of Surgery*. 30 de enero de 2013.
- Bollero D y cols. "Negative pressure surgical management after pathological scar surgical excision: a first report". *International Wound Journal*. 19 de febrero de 2013 [publicación electrónica antes de impresión]
- Bonds AM y cols. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery". *Diseases of the Colon & Rectum*. 1 de diciembre de 2013; 56(12):1403-8.
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients". *Journal of Cardiothoracic Surgery*. 6 de diciembre de 2011; 6(1):160.
- Conde-Green A y cols. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study". *Annals of Plastic Surgery*, Octubre de 2013 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown". *Journal of Wound Care*. Agosto de 2012; 21(8):386-8.
- Faroqi L y cols. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively". *Journal of Diabetic Foot Complications*. 1 de septiembre de 2013; 1;5(2):44-7.
- Glaser DA y cols. "Negative pressure therapy for closed spine incisions: A pilot study". *Wounds*. Noviembre de 2012; 24(11):308-11.
- Gomoll AH y cols. "Incisional vacuum-assisted closure therapy". *Journal of Orthopaedic Trauma*. Nov-dic de 2006; 20(10):705-9.
- Grauhan O y cols. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy". *Journal of Thoracic and Cardiovascular Surgery*. Mayo de 2013; 145(5):1387-92.
- Haghshenasskashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass". *The British Journal of Diabetes and Vascular Disease*. Ene-feb de 2011; 11(1):21-4.
- Ingariola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review". *Eplasty*. 20 de septiembre de 2013; 13:e49.
- Karlakki S y cols. "pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication". *Bone and Joint Research*. 1 de diciembre de 2013; 2(12):276-84.
- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system". *Wound Repair and Regeneration*. Septiembre de 2011;19(5):588-96.
- Mark KS y cols. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study". *Surgical Innovations* 20 de septiembre de 2013.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study". *Annals of Surgery*. Junio de 2012; 255(6):1043-7.

- Matatov T y cols. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients". Journal of Vascular Surgery. Marzo de 2013; 57(3):791-5.
- Pachowsky M y cols. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty". International Orthopaedics. Abril de 2012; 36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions". Surgical Infections. Junio de 2013; 14(3):270-4.
- Reddix RN y cols. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery". Journal of Surgical Orthopaedic Advances. Junio de 2010; 19(2):91-7.
- Stannard JP y cols. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma". Journal of Trauma. Junio de 2006; 60(6):1301-6.
- Stannard JP y cols. "Incisional negative pressure wound therapy after high-risk lower extremity fractures". Journal of Orthopedic Trauma. Junio de 2012; 26(1):37-42.
- Tauber R y cols. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications". Journal of Plastic, Reconstructive, and Aesthetic Surgery. Marzo de 2013; 66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures". American Journal of Surgery. Diciembre de 2012; 204(6):1021-4.
- Wilkes RP y cols. "Closed incision management with negative pressure wound therapy (CIM): biomechanics". Surgical Innovations. Marzo de 2012; 19(1):67-75.

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