

# Sterile Resterilizable\* Gel Breast Implant Sizers



• Resterilizable Gel Breast Implant Sizers for Volume Determination

• Shipped Sterile – \*Resterilizable up to 10 Times

• Enhanced Autoclaving Protocols

• Sizer Clearly Marked for Simple Indication of the Volume Size and the Style

• For Temporary Intraoperative Use

• Dimensions Compatible with our [MemoryGel®](#) Product Lines



Not for distribution in the U.S.



# Sterile Resterilizable\* Gel Breast Implant Sizers

## CPG™ 331<sup>NEW</sup> Cohesive III™, Tall Height, Moderate Projection

Vol.	Width	Height	Proj.	Catalog #
125 cc	9,0 cm	9,2 cm	3,2 cm	<b>RSZ-0903</b>
150 cc	9,5 cm	9,7 cm	3,3 cm	<b>RSZ-0953</b>
175 cc	10,0 cm	10,2 cm	3,4 cm	<b>RSZ-1003</b>
200 cc	10,5 cm	10,7 cm	3,5 cm	<b>RSZ-1053</b>
230 cc	11,0 cm	11,3 cm	3,6 cm	<b>RSZ-1103</b>
265 cc	11,5 cm	11,8 cm	3,7 cm	<b>RSZ-1153</b>
300 cc	12,0 cm	12,3 cm	3,9 cm	<b>RSZ-1203</b>
340 cc	12,5 cm	12,8 cm	4,0 cm	<b>RSZ-1253</b>
380 cc	13,0 cm	13,3 cm	4,1 cm	<b>RSZ-1303</b>
425 cc	13,5 cm	13,8 cm	4,3 cm	<b>RSZ-1353</b>
475 cc	14,0 cm	14,3 cm	4,5 cm	<b>RSZ-1403</b>
530 cc	14,5 cm	14,8 cm	4,6 cm	<b>RSZ-1453</b>
585 cc	15,0 cm	15,3 cm	4,8 cm	<b>RSZ-1503</b>
645 cc	15,5 cm	15,9 cm	5,0 cm	<b>RSZ-1553</b>

## CPG™ 332, Cohesive III™, Tall Height, Moderate *Plus* Projection

Vol.	Width	Height	Proj.	Catalog #
145 cc	9,0 cm	9,4 cm	3,8 cm	<b>RSZ-0909</b>
175 cc	9,5 cm	9,9 cm	4,0 cm	<b>RSZ-0959</b>
205 cc	10,0 cm	10,4 cm	4,2 cm	<b>RSZ-1009</b>
235 cc	10,5 cm	10,9 cm	4,4 cm	<b>RSZ-1059</b>
270 cc	11,0 cm	11,5 cm	4,7 cm	<b>RSZ-1109</b>
305 cc	11,5 cm	12,0 cm	4,9 cm	<b>RSZ-1159</b>
350 cc	12,0 cm	12,5 cm	5,1 cm	<b>RSZ-1209</b>
395 cc	12,5 cm	13,0 cm	5,3 cm	<b>RSZ-1259</b>
445 cc	13,0 cm	13,5 cm	5,5 cm	<b>RSZ-1309</b>
495 cc	13,5 cm	14,1 cm	5,7 cm	<b>RSZ-1359</b>
555 cc	14,0 cm	14,6 cm	5,9 cm	<b>RSZ-1409</b>
615 cc	14,5 cm	15,1 cm	6,1 cm	<b>RSZ-1459</b>
680 cc	15,0 cm	15,6 cm	6,3 cm	<b>RSZ-1509</b>

## CPG™ 333, Cohesive III™, Tall Height, High Projection

Vol.	Width	Height	Proj.	Catalog #
180 cc	9,0 cm	9,4 cm	4,6 cm	<b>RSZ-4090</b>
215 cc	9,5 cm	9,9 cm	4,8 cm	<b>RSZ-4095</b>
250 cc	10,0 cm	10,4 cm	5,1 cm	<b>RSZ-4100</b>
290 cc	10,5 cm	10,9 cm	5,3 cm	<b>RSZ-4105</b>
330 cc	11,0 cm	11,5 cm	5,6 cm	<b>RSZ-4110</b>
380 cc	11,5 cm	12,0 cm	5,8 cm	<b>RSZ-4115</b>
430 cc	12,0 cm	12,5 cm	6,0 cm	<b>RSZ-4120</b>
485 cc	12,5 cm	13,0 cm	6,2 cm	<b>RSZ-4125</b>
545 cc	13,0 cm	13,5 cm	6,5 cm	<b>RSZ-4130</b>
610 cc	13,5 cm	14,1 cm	6,7 cm	<b>RSZ-4135</b>
680 cc	14,0 cm	14,6 cm	6,9 cm	<b>RSZ-4140</b>
755 cc	14,5 cm	15,1 cm	7,1 cm	<b>RSZ-4145</b>

## Round Low Profile

Vol.	Diam.	Proj.	Catalog #
100 cc	9,4 cm	1,9 cm	<b>RSZ-1100</b>
125 cc	10,1 cm	2,0 cm	<b>RSZ-1125</b>
150 cc	10,7 cm	2,2 cm	<b>RSZ-1150</b>
175 cc	11,3 cm	2,3 cm	<b>RSZ-1175</b>
200 cc	11,8 cm	2,4 cm	<b>RSZ-1200</b>
225 cc	12,3 cm	2,5 cm	<b>RSZ-1225</b>
250 cc	12,7 cm	2,6 cm	<b>RSZ-1250</b>
275 cc	13,1 cm	2,7 cm	<b>RSZ-1275</b>
300 cc	13,5 cm	2,8 cm	<b>RSZ-1300</b>
325 cc	13,9 cm	2,9 cm	<b>RSZ-1325</b>
350 cc	14,2 cm	2,9 cm	<b>RSZ-1350</b>
375 cc	14,5 cm	2,9 cm	<b>RSZ-1375</b>
400 cc	14,8 cm	3,0 cm	<b>RSZ-1400</b>
450 cc	15,4 cm	3,1 cm	<b>RSZ-1450</b>
500 cc	16,0 cm	3,3 cm	<b>RSZ-1500</b>
550 cc	16,5 cm	3,4 cm	<b>RSZ-1550</b>
600 cc	17,0 cm	3,5 cm	<b>RSZ-1600</b>

## Round Moderate Profile

Vol.	Diam.	Proj.	Catalog #
100 cc	9,3 cm	2,1 cm	<b>RSZ-7100</b>
125 cc	10,0 cm	2,2 cm	<b>RSZ-7125</b>
150 cc	10,3 cm	2,3 cm	<b>RSZ-7150</b>
175 cc	11,2 cm	2,4 cm	<b>RSZ-7175</b>
200 cc	11,7 cm	2,5 cm	<b>RSZ-7200</b>
225 cc	12,2 cm	2,6 cm	<b>RSZ-7225</b>
250 cc	12,3 cm	2,8 cm	<b>RSZ-7250</b>
275 cc	13,2 cm	2,9 cm	<b>RSZ-7275</b>
300 cc	13,5 cm	3,0 cm	<b>RSZ-7300</b>
325 cc	13,9 cm	3,0 cm	<b>RSZ-7325</b>
350 cc	14,2 cm	3,1 cm	<b>RSZ-7350</b>
375 cc	14,4 cm	3,2 cm	<b>RSZ-7375</b>
400 cc	14,5 cm	3,2 cm	<b>RSZ-7400</b>
450 cc	14,9 cm	3,4 cm	<b>RSZ-7450</b>
500 cc	15,2 cm	3,6 cm	<b>RSZ-7500</b>
550 cc	15,9 cm	3,6 cm	<b>RSZ-7550</b>
600 cc	16,5 cm	3,7 cm	<b>RSZ-7600</b>
700 cc	17,4 cm	3,9 cm	<b>RSZ-7700</b>
800 cc	18,2 cm	4,1 cm	<b>RSZ-7800</b>

## CPG™ 321, Cohesive III™, Medium Height, Moderate Projection

Vol.	Width	Height	Proj.	Catalog #
120 cc	9,0 cm	8,5 cm	3,3 cm	<b>RSZ-0908</b>
135 cc	9,5 cm	8,9 cm	3,5 cm	<b>RSZ-0958</b>
155 cc	10,0 cm	9,4 cm	3,7 cm	<b>RSZ-1008</b>
180 cc	10,5 cm	9,9 cm	3,8 cm	<b>RSZ-1058</b>
215 cc	11,0 cm	10,3 cm	3,9 cm	<b>RSZ-1108</b>
245 cc	11,5 cm	10,8 cm	4,0 cm	<b>RSZ-1158</b>
280 cc	12,0 cm	11,3 cm	4,2 cm	<b>RSZ-1208</b>
315 cc	12,5 cm	11,8 cm	4,4 cm	<b>RSZ-1258</b>
355 cc	13,0 cm	12,2 cm	4,6 cm	<b>RSZ-1308</b>
395 cc	13,5 cm	12,7 cm	4,7 cm	<b>RSZ-1358</b>
440 cc	14,0 cm	13,2 cm	4,9 cm	<b>RSZ-1408</b>
480 cc	14,5 cm	13,6 cm	5,0 cm	<b>RSZ-1458</b>
530 cc	15,0 cm	14,1 cm	5,2 cm	<b>RSZ-1508</b>
640 cc	16,0 cm	15,0 cm	5,6 cm	<b>RSZ-1608</b>
775 cc	17,0 cm	16,0 cm	5,9 cm	<b>RSZ-1708</b>

## CPG™ 322, Cohesive III™, Medium Height, Moderate *Plus* Projection

Vol.	Width	Height	Proj.	Catalog #
140 cc	9,0 cm	8,5 cm	3,8 cm	<b>RSZ-0905</b>
165 cc	9,5 cm	8,9 cm	4,0 cm	<b>RSZ-0955</b>
195 cc	10,0 cm	9,4 cm	4,2 cm	<b>RSZ-1005</b>
225 cc	10,5 cm	9,9 cm	4,4 cm	<b>RSZ-1055</b>
255 cc	11,0 cm	10,3 cm	4,7 cm	<b>RSZ-1105</b>
295 cc	11,5 cm	10,8 cm	4,9 cm	<b>RSZ-1155</b>
330 cc	12,0 cm	11,3 cm	5,1 cm	<b>RSZ-1205</b>
375 cc	12,5 cm	11,8 cm	5,3 cm	<b>RSZ-1255</b>
420 cc	13,0 cm	12,2 cm	5,5 cm	<b>RSZ-1305</b>
475 cc	13,5 cm	12,7 cm	5,7 cm	<b>RSZ-1355</b>
525 cc	14,0 cm	13,2 cm	5,9 cm	<b>RSZ-1405</b>
585 cc	14,5 cm	13,6 cm	6,1 cm	<b>RSZ-1455</b>
650 cc	15,0 cm	14,1 cm	6,3 cm	<b>RSZ-1505</b>

## CPG™ 323, Cohesive III™, Medium Height, High Projection

Vol.	Width	Height	Proj.	Catalog #
165 cc	9,0 cm	8,5 cm	4,6 cm	<b>RSZ-0902</b>
195 cc	9,5 cm	8,9 cm	4,8 cm	<b>RSZ-0952</b>
225 cc	10,0 cm	9,4 cm	5,1 cm	<b>RSZ-1002</b>
260 cc	10,5 cm	9,9 cm	5,3 cm	<b>RSZ-1052</b>
300 cc	11,0 cm	10,3 cm	5,6 cm	<b>RSZ-1102</b>
345 cc	11,5 cm	10,8 cm	5,8 cm	<b>RSZ-1152</b>
390 cc	12,0 cm	11,3 cm	6,0 cm	<b>RSZ-1202</b>
440 cc	12,5 cm	11,8 cm	6,2 cm	<b>RSZ-1252</b>
495 cc	13,0 cm	12,2 cm	6,5 cm	<b>RSZ-1302</b>
555 cc	13,5 cm	12,7 cm	6,7 cm	<b>RSZ-1352</b>
620 cc	14,0 cm	13,2 cm	6,9 cm	<b>RSZ-1402</b>
685 cc	14,5 cm	13,6 cm	7,1 cm	<b>RSZ-1452</b>

## Round Moderate *Plus* Profile

Vol.	Diam.	Proj.	Catalog #
100 cc	8,2 cm	2,7 cm	<b>RSZ-1001</b>
125 cc	8,9 cm	2,8 cm	<b>RSZ-1251</b>
150 cc	9,5 cm	2,9 cm	<b>RSZ-1501</b>
175 cc	10,0 cm	3,1 cm	<b>RSZ-1751</b>
200 cc	10,5 cm	3,2 cm	<b>RSZ-2001</b>
225 cc	10,9 cm	3,3 cm	<b>RSZ-2251</b>
250 cc	11,3 cm	3,4 cm	<b>RSZ-2501</b>
275 cc	11,7 cm	3,5 cm	<b>RSZ-2751</b>
300 cc	12,0 cm	3,6 cm	<b>RSZ-3001</b>
325 cc	12,3 cm	3,8 cm	<b>RSZ-3251</b>
350 cc	12,5 cm	3,9 cm	<b>RSZ-3501</b>
375 cc	12,8 cm	4,0 cm	<b>RSZ-3751</b>
400 cc	13,1 cm	4,0 cm	<b>RSZ-4001</b>
450 cc	13,6 cm	4,2 cm	<b>RSZ-4501</b>
500 cc	14,1 cm	4,3 cm	<b>RSZ-5001</b>
550 cc	14,6 cm	4,5 cm	<b>RSZ-5501</b>
600 cc	15,0 cm	4,6 cm	<b>RSZ-6001</b>
700 cc	15,8 cm	4,9 cm	<b>RSZ-7001</b>
800 cc	16,5 cm	5,1 cm	<b>RSZ-8001</b>

## Round High Profile

Vol.	Diam.	Proj.	Catalog #
125 cc	8,3 cm	3,5 cm	<b>RSZ-1254</b>
150 cc	8,8 cm	3,7 cm	<b>RSZ-1504</b>
175 cc	9,3 cm	3,9 cm	<b>RSZ-1754</b>
200 cc	9,7 cm	4,0 cm	<b>RSZ-2004</b>
225 cc	10,1 cm	4,2 cm	<b>RSZ-2254</b>
250 cc	10,5 cm	4,3 cm	<b>RSZ-2504</b>
275 cc	10,8 cm	4,4 cm	<b>RSZ-2754</b>
300 cc	11,1 cm	4,5 cm	<b>RSZ-3004</b>
325 cc	11,4 cm	4,6 cm	<b>RSZ-3254</b>
350 cc	11,7 cm	4,8 cm	<b>RSZ-3504</b>
375 cc	12,0 cm	4,8 cm	<b>RSZ-3754</b>
400 cc	12,2 cm	5,0 cm	<b>RSZ-4004</b>
425 cc	12,5 cm	5,0 cm	<b>RSZ-4254</b>
450 cc	12,8 cm	5,1 cm	<b>RSZ-4504</b>
500 cc	13,2 cm	5,3 cm	<b>RSZ-5004</b>
550 cc	13,6 cm	5,5 cm	<b>RSZ-5504</b>
600 cc	14,0 cm	5,6 cm	<b>RSZ-6004</b>
650 cc	14,4 cm	5,7 cm	<b>RSZ-6504</b>
700 cc	14,8 cm	5,8 cm	<b>RSZ-7004</b>
800 cc	15,5 cm	6,0 cm	<b>RSZ-8004</b>

## CPG™ 311<sup>NEW</sup> Cohesive III™, Low Height, Moderate Projection

Vol.	Width	Height	Proj.	Catalog #
100 cc	9,0 cm	8,0 cm	3,0 cm	<b>RSZ-0090</b>
120 cc	9,5 cm	8,5 cm	3,2 cm	<b>RSZ-0095</b>
140 cc	10,0 cm	8,9 cm	3,4 cm	<b>RSZ-0100</b>
160 cc	10,5 cm	9,4 cm	3,6 cm	<b>RSZ-0105</b>
180 cc	11,0 cm	9,8 cm	3,7 cm	<b>RSZ-0110</b>
210 cc	11,5 cm	10,3 cm	3,9 cm	<b>RSZ-0115</b>
235 cc	12,0 cm	10,7 cm	4,1 cm	<b>RSZ-0120</b>
270 cc	12,5 cm	11,2 cm	4,2 cm	<b>RSZ-0125</b>
300 cc	13,0 cm	11,6 cm	4,4 cm	<b>RSZ-0130</b>
335 cc	13,5 cm	12,0 cm	4,6 cm	<b>RSZ-0135</b>
375 cc	14,0 cm	12,5 cm	4,7 cm	<b>RSZ-0140</b>
415 cc	14,5 cm	12,9 cm	4,9 cm	<b>RSZ-0145</b>
460 cc	15,0 cm	13,4 cm	5,1 cm	<b>RSZ-0150</b>
510 cc	15,5 cm	13,8 cm	5,2 cm	<b>RSZ-0155</b>
560 cc	16,0 cm	14,3 cm	5,4 cm	<b>RSZ-0160</b>
615 cc	16,5 cm	14,7 cm	5,6 cm	<b>RSZ-0165</b>

## CPG™ 312, Cohesive III™, Low Height, Moderate *Plus* Projection

Vol.	Width	Height	Proj.	Catalog #
125 cc	9,0 cm	8,0 cm	3,8 cm	<b>RSZ-0907</b>
145 cc	9,5 cm	8,4 cm	4,0 cm	<b>RSZ-0957</b>
170 cc	10,0 cm	8,8 cm	4,2 cm	<b>RSZ-1007</b>
195 cc	10,5 cm	9,3 cm	4,4 cm	<b>RSZ-1057</b>
225 cc	11,0 cm	9,7 cm	4,7 cm	<b>RSZ-1107</b>
255 cc	11,5 cm	10,2 cm	4,9 cm	<b>RSZ-1157</b>
290 cc	12,0 cm	10,6 cm	5,1 cm	<b>RSZ-1207</b>
330 cc	12,5 cm	11,1 cm	5,3 cm	<b>RSZ-1257</b>
370 cc	13,0 cm	11,5 cm	5,5 cm	<b>RSZ-1307</b>
415 cc	13,5 cm	12,0 cm	5,7 cm	<b>RSZ-1357</b>
465 cc	14,0 cm	12,4 cm	5,9 cm	<b>RSZ-1407</b>
515 cc	14,5 cm	12,8 cm	6,1 cm	<b>RSZ-1457</b>
570 cc	15,0 cm	13,3 cm	6,3 cm	<b>RSZ-1507</b>
690 cc	16,0 cm	14,2 cm	6,8 cm	<b>RSZ-1607</b>

## CPG™ 313<sup>NEW</sup> Cohesive III™, Low Height, High Projection

Vol.	Width	Height	Proj.	Catalog #
130 cc	9,0 cm	8,1 cm	4,4 cm	<b>RSZ-0906</b>
155 cc	9,5 cm	8,6 cm	4,5 cm	<b>RSZ-0956</b>
180 cc	10,0 cm	9,0 cm	4,6 cm	<b>RSZ-1006</b>
210 cc	10,5 cm	9,5 cm	4,8 cm	<b>RSZ-1056</b>
240 cc	11,0 cm	9,9 cm	4,9 cm	