

JIN TYPE - FEMORAL



JIN TYPE - FEMORAL NAILING SYSTEM



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & Instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and innovator products that exceed the highest expectations of our customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 60 countries worldwide.

Our Achievements

US FDA 510K APPROVED •
(Locking Plates and Screws)

EC Certificate - Class IIa & IIb •

ISO 13485:2016 •

ISO 9001:2015 •

India FDA Approved •

GMP Certificate •

Orthopaedic Implants

www.auxein.com


JIN Type -Femoral Nail

Ø Dia	Length	Stainless Steel	Titanium
9.0mm	320mm	1293-9-320	TI-1293-9-320
9.0mm	340mm	1293-9-340	TI-1293-9-340
9.0mm	360mm	1293-9-360	TI-1293-9-360
9.0mm	380mm	1293-9-380	TI-1293-9-380
9.0mm	400mm	1293-9-400	TI-1293-9-400
9.0mm	420mm	1293-9-420	TI-1293-9-420
9.0mm	440mm	1293-9-440	TI-1293-9-440
10.0mm	320mm	1293-10-320	TI-1293-10-320
10.0mm	340mm	1293-10-340	TI-1293-10-340
10.0mm	360mm	1293-10-360	TI-1293-10-360
10.0mm	380mm	1293-10-380	TI-1293-10-380
10.0mm	400mm	1293-10-400	TI-1293-10-400
10.0mm	420mm	1293-10-420	TI-1293-10-420
10.0mm	440mm	1293-10-440	TI-1293-10-440
11.0mm	320mm	1293-11-320	TI-1293-11-320
11.0mm	340mm	1293-11-340	TI-1293-11-340
11.0mm	360mm	1293-11-360	TI-1293-11-360
11.0mm	380mm	1293-11-380	TI-1293-11-380
11.0mm	400mm	1293-11-400	TI-1293-11-400
11.0mm	420mm	1293-11-420	TI-1293-11-420
11.0mm	440mm	1293-11-440	TI-1293-11-440
12.0mm	320mm	1293-12-320	TI-1293-12-320
12.0mm	340mm	1293-12-340	TI-1293-12-340
12.0mm	360mm	1293-12-360	TI-1293-12-360
12.0mm	380mm	1293-12-380	TI-1293-12-380
12.0mm	400mm	1293-12-400	TI-1293-12-400
12.0mm	420mm	1293-12-420	TI-1293-12-420
12.0mm	440mm	1293-12-440	TI-1293-12-440
13.0mm	320mm	1293-13-320	TI-1293-13-320
13.0mm	340mm	1293-13-340	TI-1293-13-340
13.0mm	360mm	1293-13-360	TI-1293-13-360
13.0mm	380mm	1293-13-380	TI-1293-13-380
13.0mm	400mm	1293-13-400	TI-1293-13-400
13.0mm	420mm	1293-13-420	TI-1293-13-420
13.0mm	440mm	1293-13-440	TI-1293-13-440

**Ø4.5mm Locking Bolt, Self Tapping, For JIN Type -Femoral Nail**

Length	Stainless Steel	Titanium
20mm	1293-4.5-20	TI-1293-4.5-20
25mm	1293-4.5-25	TI-1293-4.5-25
30mm	1293-4.5-30	TI-1293-4.5-30
35mm	1293-4.5-35	TI-1293-4.5-35
40mm	1293-4.5-40	TI-1293-4.5-40
45mm	1293-4.5-45	TI-1293-4.5-45
50mm	1293-4.5-50	TI-1293-4.5-50
55mm	1293-4.5-55	TI-1293-4.5-55
60mm	1293-4.5-60	TI-1293-4.5-60
65mm	1293-4.5-65	TI-1293-4.5-65
70mm	1293-4.5-70	TI-1293-4.5-70
75mm	1293-4.5-75	TI-1293-4.5-75
80mm	1293-4.5-80	TI-1293-4.5-80

**1293-001 End Cap for JIN Type - Femoral Nail**

INS-1293-01 Drill Sleeve 8/4mm For Drill Bit Ø4mm



INS-1293-02 Protection Sleeve 10/8.2mm



INS-1293-03 Drill Sleeve 8/5.3mm For Drill Bit Ø5.2mm



INS-1293-04 Tissue Protector



INS-1293-05 Right Angle Locking Rod Ø4mm For Nail Holding Bolt



INS-1293-06 T-Handle Screw Driver, Hex 3.5mm For Ø4.5mm Locking Bolt, Self Tapping



INS-1293-07 Positioning Flat Drill Ø5.2mm For Positioning Rod



INS-1293-08 T-Handle Tap Ø6mm For Ø4.5mm Locking Bolt, Self Tapping



INS-1293-09 Screws Exactor Ø9.9mm



INS-1293-10 Positioning Rod Ø5mm



INS-1293-12 Allen Wrench, Hex 5mm For Distal/Proximal Bolt



INS-1293-13 Allen Wrench, Hex 3mm For Drill Bit Stopper Ø4.0mm



INS-1293-14 Combination Wrench, 12mm



INS-1293-15 Cannulated Bone Awl 12/4mm



Reamer Head



Code	Dia
INS-1293-24	Ø13mm
INS-1293-30	Ø12.5mm
INS-1293-16	Ø12mm
INS-1293-17	Ø11.5mm
INS-1293-18	Ø11mm
INS-1293-19	Ø10.5mm
INS-1293-20	Ø10mm
INS-1293-21	Ø9.5mm
INS-1293-22	Ø9mm
INS-1293-42	Ø8.5mm
INS-1293-43	Ø8mm

INS-1293-23 Drill Bit Plain Shank/ Jacob Chuck End Ø4mm x Length 300mm
For Ø4.5mm Locking Bolt, Self Tapping



INS-1293-27 Stopper For Drill Bit Ø4mm



INS-1293-25 Drill Bit Plain Shank/ Jacob Chuck End Ø5.2mm x Length 250mm




INS-1293-26 Entry Reamer Ø12.5mm With Quick Coupling End



INS-1293-28 Reducer, Cannulation 4.3mm



INS-1293-29 Flexible Reaming Shaft (For Reamer Head Ø8mm to 13mm) (For 2.5mm Wire Pass)



INS-1293-31 Insertion Handle For JIN Type -Femoral Nail



INS-1293-32

Distal Targeting Arm

**INS-1293-33**

T-Handle Quick Coupling For Flexible Shaft/Entry Reamer

**INS-1293-34**

Positioning Block

**INS-1293-35**

Trocar Ø8mm

**INS-1293-36**

Distal Targeting Device (Length 320mm To 420mm)

**INS-1293-30**

Bolt For Proximal Targetting Device



INS-1293-37

Drill Sleeve 15.5/13.5mm For Entry Reamer

**INS-1293-38**

Measuring Device upto 90mm For Ø4.5mm Locking Bolt, Self Tapping

**INS-1293-39**

Extraction Rod With Ram

**INS-1293-40**

Flexible Universal Adapter Ø7.9mm

**INS-1293-41**

Impactor Ø9.9mm

**INS-1293-44**

Guide Wire with Ball Tip Ø3.5mm/2.5mm x Length 1050mm - Nitinol



INS-1293-46

Nail Holding Bolt

**INS-1293-47**

Bolt For Proximal Targeting Device

**INS-1293-48**

Proximal Targeting Device

**INS-1293-49**Guide Wire Introducer
(For 2.5mm Wire)**7-022-01**

Radiographic Ruler For Long JIN-Type Femoral Nail

**7-022-02**

Plastic Exchange Tube For JIN-Type Femoral Nail



INS-1293T



Trays For JIN-Type Femoral Nail

**INS-1293B**

Container With Lid



INS-1293 Instrument Set For JIN- Type Femoral Nail





INS-1293 Instrument Set For JIN-Type Femoral Nail

Code	Set Consisting of	Units
INS-1293-01	Drill Sleeve 8/4mm For Drill Bit Ø4mm	2
INS-1293-02	Protection Sleeve 10/8.2mm	3
INS-1293-03	Drill Sleeve 8/5.3mm For Drill Bit Ø5.2mm	1
INS-1293-04	Tissue protector	1
INS-1293-05	Right Angle Locking Rod Ø4mm For Nail Holding Bolt	1
INS-1293-06	T-Handle Screw Driver, Hex 3.5mm For Ø4.5mm Locking Bolt, Self Tapping	1
INS-1293-07	Positioning Flat Drill Ø5.2mm For Positioning Rod	1
INS-1293-08	T-Handle Tap Ø6mm For Ø4.5mm Locking Bolt, Self Tapping	1
INS-1293-09	Screws Exactor Ø9.9mm	1
INS-1293-10	Positioning Rod Ø5mm	1
INS-1293-12	Allen Wrench, Hex 5mm For Distal/Proximal Bolt	1
INS-1293-13	Allen Wrench, Hex 3mm For Drill Bit Stopper Ø4.0mm	1
INS-1293-14	Combination Wrench, 12mm	1
INS-1293-15	Cannulated Bone Awl 12/4mm	1
INS-1293-24	Reamer Head Ø13mm	1
INS-1293-30	Reamer Head Ø12.5mm	1
INS-1293-16	Reamer Head Ø12mm	1
INS-1293-17	Reamer Head Ø11.5mm	1
INS-1293-18	Reamer Head Ø11mm	1
INS-1293-19	Reamer Head Ø10.5mm	1
INS-1293-20	Reamer Head Ø10mm	1
INS-1293-21	Reamer Head Ø9.5mm	1
INS-1293-22	Reamer Head Ø9mm	1
INS-1293-42	Reamer Head Ø8.5mm	1
INS-1293-43	Reamer Head Ø8mm	1
INS-1293-23	Drill Bit Plain Shank/ Jacob Chuck End Ø4mm x Length 300mm For Ø4.5mm Locking Bolt, Self Tapping	2
INS-1293-27	Stopper For Drill Bit Ø4mm	1
INS-1293-25	Drill Bit Plain Shank/ Jacob Chuck End Ø5.2mm x Length 250mm	1
INS-1293-26	Entry Reamer Ø12.5mm With Quick Coupling End	1
INS-1293-28	Reducer, Cannulation 4.3mm	1
INS-1293-29	Flexible Reaming Shaft (For Reamer Head Ø8mm to 13mm) (For 2.5mm Wire Pass)	2
INS-1293-31	Insertion Handle For JIN Type -Femoral Nail	1
INS-1293-32	Distal Targeting Arm	1
INS-1293-33	T-Handle Quick Coupling For Flexible Shaft/Entry Reamer	1

Code	Set Consisting of	Units
INS-1293-34	Positioning Block	1
INS-1293-35	Trocars Ø8mm	1
INS-1293-36	Distal Targeting Device (Length 320mm To 420mm)	1
INS-1293-30	Bolt For Distal Targeting Device	3
INS-1293-37	Drill Sleeve 15.5/13.5mm For Entry Reamer	1
INS-1293-38	Measuring Device upto 90mm For Ø4.5mm Locking Bolt, Self Tapping	1
INS-1293-39	Extraction Rod With Ram	1
INS-1293-40	Flexible Universal Adapter Ø7.9mm	1
INS-1293-41	Impactor Ø9.9mm	1
INS-1293-44	Guide Wire with Ball Tip Ø3.5mm/2.5mm x Length 1050mm - Nitinol	2
INS-1293-46	Nail Holding Bolt	1
INS-1293-47	Bolt For Proximal Targeting Device	1
INS-1293-48	Proximal Targeting Device	1
INS-1293-49	Guide Wire Introducer(For 2.5mm Wire)	1
7-022-01	Radiographic Ruler For Long JIN-Type Femoral Nail	1
7-022-02	Plastic Exchange Tube For JIN-Type Femoral Nail	1
INS-1293T	Trays For JIN- Type Femoral Nail	1
INS-1293B	Container With Lid	1

7-006

Flexible Reamer Instrument set


1459-006

Flexible Reamer Shaft (AO Connection)


Reamer Head

Code	Dia
7-006-01	Ø8mm
1459-008	Ø8.5mm
1459-009	Ø9mm
1459-010	Ø9.5mm
1459-011	Ø10mm
1459-012	Ø10.5mm
1459-013	Ø11mm
1459-014	Ø11.5mm

Code	Dia
1459-015	Ø12mm
1459-016	Ø12.5mm
1459-017	Ø13mm
1459-018	Ø13.5mm
1459-019	Ø14mm
7-006-02	Ø14.5mm
7-006-03	Ø15mm



**Corporate Office**

7/A, Nirankari Colony, Delhi-110009, India

Manufacturing

Plot No. 169-170, Phase-IV, Sector-57, Kundli Industrial Area,
HSIIDC, Sonepat-131028, Haryana, India

Tel: +91 99106 43638 | Fax: +91 86077 70197

Email: info@auxein.com | Web: www.auxein.com